


HOW TO RETURN SELF-COLLECTED SWABS AND URINE SAMPLES BY MAIL


STEP 1

- Place your completed samples into their own resealable biohazard plastic bags (labelled Bag #1).
- Pinch the zipper on each Bag #1 to seal closed.




STEP 2

- Place sealed Bag #1 containing the sample tube into the plastic bag labelled Bag #2.
- If you are returning more than one sample, put all Bag #1s into the same Bag #2.
- Do not remove the absorbent pack from Bag #2.



STEP 3





Seal Bag #2:

1. Peel off the adhesive strip at the top of Bag #2.
2. Fold the top of the bag downward on the dotted line to seal the bag.

STEP 4

- Place sealed Bag #2 into the shipping box.
- Close the box lid.



STEP 5



Adhesive

- Seal the shipping box by placing the adhesive strip over the opening of the box.

STEP 6

- Place the shipping box into a post office box (mail receptacle) or bring it to a local post office.
- You will not be charged shipping for mailing this box.

